pedalo®

Rolli-Parcours Wheelchair Course

20250000

www.pedalo.de

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Driving a wheelchair is not easy!

In everyday life, wheelchair users still deal with obstacles such as curbsides, short steep inclines, getting on and off bus stops, stairways, inclined planes or entrances to shopping centers for example. This requires a lot of skill but also the courage to even try. Inexperienced wheelchair users often quickly reach their personal limits and find themselves in dangerous situations.

Thus, Pedalo has developed a wheelchair course for technique training together with wheelchair users to help them find their way around such everyday situations.

The Wheelchair Course

Together with wheelchair users Pedalo has developed a wheelchair course for technique training which helps finding a way around such everyday situations. The course enables them to adapt and improve their skills and behavior to the relevant situations. At the same time, new movement experiences are learned, such as rocking by using the upper body, which also helps them to get better skills.

Great benefits

- Technique training improves the handling of the wheelchair
- Balance training provides more safety on uneven surfaces
- Learning to take steps, climbs, descents and inclined planes better

- Fun and motivation
- Can also be used for road safety education and safety training for cyclists and inline skaters

The advantages in use

- Quick, tool-free assembly and disassembly (plug-in system)
- Can be handled by one person
- Levels of difficulty can be individually adjusted
- Mobile use, fits in cars with a loading area of 180 cm in length

The course consists of 5 stations

Seesaw – Inclined Plane – Vibration Plate – Turning Plate – Step Plate

In use at the German Red Cross KV Ulm e.V. in cooperation with Habila

The everyday life of people in wheelchairs is full of challenges that are often unclear to people with no experience of wheelchairs or the subject area. The wheelchair course is a very good way to illustrate precisely this aspect. The project, which aims to give interested people an insight into the hurdles and possibilities of wheelchair use, offers the option of getting to know exactly this, experiencing it independently and trying it out.

Wheelchair Course Station Inclined Plane


Roads, sidewalks and properties are not always flat and level. If wheelchair users get into a sloping position with their wheelchair, there is a risk of tipping over. This must be counteracted by the upper body. When driving on the inclined plane (slope), the wheelchair's driving behavior is changed in such a way that it very quickly steers into the fall line – if the ride is uneven. By using this station, wheelchair users learn to master also this situation very quickly.

What makes this station with its inclined plane such special

Rubberized up and down wedges provide safe access to the inclined plane. Depending on the choice of the pre-drawn ride lines, riders can take different inclined positions.

Specification: Water-repellent, black-coated sloping boards with non-slip surface finish, can be assembled and disassembled without tools. Consisting of 3 birch veneer panels 148x49x2.1cm and 3 insert boards 16cm high. Plus 2 insertable edge limiters 94.5x5cm and 2 rubberized grippy up and down wedges 100x19cm. The 3 inclined planes are connected by 4 tensioning straps. Complete inclined plane assembled 450x62.5x20cm. Transport dimensions per wooden board approx. 150x62.5x20cm (150x50x20cm when dismantled). Weight capacity up to 150kg. 51kg.

Exercise: Drive at an angle and try not to leave the line.

Wheelchair Course Station Seesaw

Driving up a ramp requires increased effort. A ramp that tips over halfway require an additional sense of balance. The driver and wheelchair must work in harmony to achieve balance with the wheelchair and balance the seesaw. Riding down and leaving the seesaw also requires good skill.

Specification: Water-repellent, black-coated seesaw with non-slip surface finish. Consisting of 2 birch veneer panels 148x49x2.1cm and 2 insertable edge limiters 94.5x5cm. Plus 2 pairs of anti-slip runners K70 with soft edges, 7cm high and 4 rubber buffers for low-noise contact with the floor. 2-piece rubberized grippy ramp wedge (49x19cm each). The two wooden boards are connected in parallel using 2 tensioning straps. Seesaw assembled 162x98x12cm. Can be easily dismantled into 2 individual seesaws (each 162x49x12cm) for transportation. Weight capacity up to 120kg. 27.2kg.

Exercise: Drive onto the seesaw and try to keep your balance, then carefully ride off the seesaw. The same exercise can also be done in backwards.

Wheelchair Course Station Vibration Plate

Potholes in the ground, objects on the road, rough terrain, gravel paths, unpaved cycle paths or grassy areas are often unmanageable tasks for wheelchair users. They make driving extremely difficult, and some wheelchair users feel like being on a bumpy track and are overwhelmed by the demands. It is often the many, small and nerve-wracking obstacles that annoy and often make one almost despair.

Variability of the insertable wooden strips
On the vibrating plate, there are several of these obstacles in the form of insertable wooden strips.

The stations of the wheelchair course

Up to seven obstacles per wheel, each up to 2 cm high, can make the ride very bumpy. It can even happen that almost all the wheels block at the same time. Wheelchair users learn how to free themselves from such situations by acting purposefully. The two boards can be laid parallel next to each other or shifted next to each other. With up and down wedge.

Specification: Water-repellent, black-coated vibrating plate with non-slip surface finish and rubberized anti-slip feet. Consisting of 2 birch veneer panels 148x49x2.1cm, 14 insert strips 41x4cm and 4 rubberized grippy up and down wedges 49x19cm. The two boards are simply laid next to each other without a connection. Vibration plate completely assembled 176x98x8cm (if placed parallel). Transport dimensions per wooden board with wedges but without insertion strips 176x49x5.5cm. Weight capacity up to 150kg. 30.5kg.

Exercise: Drive over the vibratory plate and pass over the obstacles (the inserted strips).

Wheelchair Course Station Turning Plate

Safe turning and maneuvering in confined spaces without bumping into things is an important requirement, especially in stores. The wheelchair turning plate is an ideal training tool for this. Riding up over the ramp wedge, turning around the own body axis in a three-sided confined space without touching the adjustable edge limiters is an interesting challenge.

Specification: Water-repellent, black-coated turning plate with non-slip surface finish and rubberized anti-slip feet. Consisting of 2 birch veneer panels 106x50x2.1cm, 3 edge limiters 94.5x5cm (can be inserted laterally in variable positions) and 1 rubberized grippy up and down wedge 100x19cm. The two wooden boards are connected using 2 tensioning straps. Turning plate assembled 113.5x106x6.5cm. Transport dimensions per board without edge strips approx. 50(64)106x3cm. Weight capacity up to 150kg. 19.2kg.

Exercise: Drive over the ramp wedge onto the turning plate. Turn around your own axis, trying not to touch the edge limiters. Then drive back down over the wedge.

The stations of the wheelchair course

Wheelchair Course Station Step Plate

Even small bumps in the floor of just a few centimeters are a real obstacle for the front wheels of a wheelchair. This is a problem that many wheelchair users have to deal with on a daily basis. These bumps or steps can only be crossed safely with a great deal of skill. The step plate can be passed and driven over from both sides. Two steps of different heights thus offer two levels of difficulty.

Specification: Water-repellent, black-coated step plate with non-slip surface finish. Anti-slip feet on the underside make it difficult to the plate slipping away. Consisting of 1 birch veneer panel 98x62x2.1 cm and 8 rubberized anti-slip feet in 2 different heights for attaching to the underside of the board. Step plate assembled 98x98x5.6(3.6)cm. Weight capacity up to 150kg. 9kg.

Exercise: Drive over the step plate from both sides and practice mastering the steps with increasing skill.

All stations are also available separately
Seesaw (20251000) – Inclined Plane (20252000) – Vibration Plate (20253000) –
Turning Plate (20254000) – Step Plate (20255000)

Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore signs of using original wood and no quality reduction.

Please note:

WARNING! The devices are only to be used for the purposes described in this user manual.

WARNING! Place the unit only on flat level ground.

After longer or frequent use, check that the screws and all installed parts are tight and retighten if necessary. Defective and worn parts should be replaced by original replacement parts immediately.

Protect item for damp and wet. This device is meant for using from persons at the age of 14 up. Maximum weight capacity approx. 150 kg

Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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