
# pedalo®


**Fußwerkstatt**  
**Foot Workshop**

16002000

[www.pedalo.de](http://www.pedalo.de)

<b>Ossa digiti pedis</b> Zehenknochen Phalanges	<b>Phalanges proximales</b> Zehengrundglieder Proximal phalanges
<b>Phalanges distales</b> Zehenendglieder Distal phalanges	<b>Metatarsalia</b> Mittelfußknochen Metatarsal bone
<b>Phalanges mediae</b> Zehennittelglieder Middle phalanges	<b>Ossa tarsi</b> Fußwurzelknochen Tarsal bone

<b>1. Talus</b> Sprungbein Ankle bone	<b>4. Os cuneiforme I</b> Keilbein Sphenoid bone	
<b>2. Calcaneus</b> Fersenbein Heel bone	<b>5. Os cuneiforme II</b> Keilbein Sphenoid bone	
<b>3. Os naviculare</b> Kahnbein Scaphoid bone	<b>6. Os cuneiforme III</b> Keilbein Sphenoid bone	<b>7. Os cuboideum</b> Würfelbein Cuboid bone

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
The human foot is a fantastic construction of 26 bones, from many ligaments, tendons and nerves. If healthy, it carries (as a pair) our full body weight, painlessly and cushions each step perfectly. Our feet regulate our balance and provide us a wide range of movement possibilities. Only with healthy feet, we have the necessary and important security of movement from a young age, which we need throughout our life for fall prevention – reactively but also preventively.

Much of the population is suffering from foot problems of which two-thirds are often caused by (wrong) movement habits and inappropriate biomechanical stress on muscles and joints, what leads to skew-, flat- or splayfeet. This inevitably means an overstressing of foot, knee and hip joints and spine. Often this is the beginning of a chain of larger foot problems and in its further course it can lead to chronic inflammation and increased signs of wear and tear in the joints, associated with pain. Shoe inlays, prescribed by a physician, support the foot position and may relieve pain but they also place the foot in its function passively and do not contribute to strengthen the insufficient foot musculature.

We at Pedalo therefore recommend a daily dose of training for your feet. In our opinion this is the only way to strengthen the arch of the foot and even to improve incipient foot deformities. Foot training begins with bare-foot walking. Then it includes a daily mobilization of feet, associated with an active muscular structure workout.

Only with healthy feet, we have fun in the movement and are able to stand, walk, run, jump, dance and bounce. To achieve all these tasks, the foot relies on a healthy and resilient vault structure.

- S9 *Plantar Stretcher*
- S1 *Twister*
- S2 *Toe Grip*
- S3 *Senso*
- S4 *Seesaw*
- S5 *Top*
- S7 *Forefoot Lifter*
- S6 *Heel Lifter*
- S8 *Foot Arch Stretcher*

In everyday life and in sports, healthy feet are a blessing and the basis for pain-free movement. Almost all malformations of feet are due to muscular deficiencies or imbalance of the foot and leg muscles—and this often already from an early age on. To strengthen the muscles a continuous training is essential. The nine exercise stations of the Pedalo-Foot Workshop are convincing through their versatility and activate all systems and structures, which are important for a healthy and functional foot. Thus, the different stations are intended to activate the feet, sensitize and massage the foot soles, stabilize the ankle joints, strengthen the foot and lower leg musculature as well as to stretch the foot arch and the plantar fascia. The Pedalo-Foot Workshop is suited ideally as training circuit for prevention and rehabilitation. It is used in fitness and health studios, podiatry, physiotherapy, and rehabilitation facilities; in retirement homes or wellness hotels; and also as part of corporate health management. The storage stand offers a space-saving and clear possibility for storage.

**Specification:** 1 storage stand 47x40x110cm, solid beech, base plate CDF, colored black. 9 exercise stations made of CDF and solid beech or birch plywood. Each integrated on base plate CDF 47x40cm, with non-slip feet. Incl. exercise pictures for the boards. The height of the stations varies between 1.2 and 17 cm depending on the item. 34kg. **16002000**

## S1 Twister 16002010

When doing the Twist, we perform a dynamic movement, and almost all muscles of the foot and the lower leg are involved. As a starting position both feet are placed side by side. Now the heel is lifted, turned on the forefoot to side and lowered down again. Then the forefoot is lifted, turned on the heel side and put down again. Repeat the same exercise in opposite direction now. Perform the exercise in a slow and controlled manner.

**Exercise scope:** Repeat 10 times left and right.

## S2 Toe Grip 16002020

Toe gripping strengthens and relaxes the foot and the toe muscles. Similar to the movement of a caterpillar, the feet move forward by stretching and contracting the toes. The cross grooves support the forward movement. Make sure that all toes grab in the grooves. Then carry out the contrary toe movement. To this extend the toes as far and wide as possible upward after each small forward movement. Perform the exercise slowly and evenly so that you can feel and move at best each single toe.

**Exercise scope:** Repeat 5 times.

## S3 Senso 16002030

The Senso board sensitizes and massages the soles of the feet, acting like acupressure. In this way, the senso station has a positive influence on the nervous system and helps to activate the metabolism. The slow and steady exercises stimulate and relax the foot muscles at the same time.

**Exercise scope:** Stand on the wooden nubs for 3 minutes but change the foot position several times, so that the nubs can constantly provide new stimuli and stimulate blood circulation.

### S4 Seesaw 16002040

Footboards Seesaw offer a linear movement of the foot into pronation and supination (strain of inner/outer side of the foot). Beside of the foot muscles, the stabilizing musculature of the lower leg is activated. The combination of these muscles influences our posture stability.

**Exercise scope:** Standing with your feet parallel, rock side-to-side several times, then stabilize the boards horizontally after a few seconds.

### S5 Top 16002050

Footboards Top tilt on all sides and are very close to the function of the ankle joint and its mobility. The requirements to the coordination skills are here higher than with the Footboards Seesaw. Therefore, this station serves for mobilization and the stabilization of the ankle joint.

**Exercise scope:** Train the ankle joint by doing several alternating circular movements left and right (clockwise, counterclockwise or both at the same time) against the resistance of the Footboards Top. In between, stop the movement for a few seconds and stabilize in a horizontal position if possible (aim for about 10 seconds). Or, as an exercise variation, stand alternately on one leg right/left.

### S6 Heel Lifter 16002060

Exercises on the Heel lifter activate the interaction of foot muscles with the lower leg muscles and give the ankle joint its necessary stability. The starting position alone (forefoot placed on the round wood, heel on the gray plate) provides a stretch in the Achilles tendon/calf muscles. By lifting the heel to a horizontal position, the muscles around the foot, ankle, and lower leg—as well as the sense of balance and responsiveness—are particularly challenged.

**Exercise scope:** Repeat 10 times.

### S7 Forefoot Lifter 16002070

People with foot weakness must have much concentration not to fail in taking steps. Every little bump means a danger zone. So now, the normally almost automated movement process of walking, represents for the person concerned—in the truest sense of the word—a task, that he must perform step by step attentively. Unfortunately this always affects his perception of the environment adversely. Very often people suffering from weak foot dorsiflexion, place their foot not with the heel first but with the forefoot or the entire sole of the foot. A natural process of rolling of the foot and forward swing of the leg is no longer possible for foot weak dorsiflexion. The person concerned becomes cramped, walks abnormally and sometimes also even with clenched toes. Swinging the leg forward is performed by a circular hip movement (Circumduction). In the long term, people who move this way complain of pain in the hip and pelvis. It is therefore necessary to counteract this as best and as quickly as possible. Anatomy: The muscle tibialis anterior supports the dorsal flexion, respectively extension, supination, and adduction of the foot. In addition, his tension stabilizes the ankle joint, especially in the first phase of ground contact when walking or running. He is the antagonist of the muscle 'peroneus longus'. The paralysis of the muscle 'tibialis anterior' leads to the so called steppage gait.

For a better gait pattern, targeted exercises can be used to activate the neuronal control of the muscles responsible for rolling the foot (mainly the anterior tibial muscles).

**Exercise:** The exercise is performed without shoes, as this is the only way to optimally engage all the structures involved in the movement. The forefoot is positioned that way so that the resistance band runs across the base joints of the toes. The toes extend beyond the resistance band, and the heel is placed on the half-round wood.

Depending on the number of resistance bands selected (2 bands are available), the resistance of the pulling force can be individually selected from level 1 to 4 (two resistance options per band).

Depending on the desired length of pull, the knee of the training foot is positioned vertically above or behind the heel.

The leg remains completely relaxed during the exercise. Concentrate on the forefoot and the front shin muscle. Make sure the toes remain extended during the exercise and are not pulled upward.

**Exercise scope:** Perform 3 sets of 10 repetitions per foot. It is important to complete all repetitions; therefore, please select the number of resistance bands accordingly so that the exercise series can be performed correctly. Do not change exercise foot during the individual exercise series. It means: First with the right foot: 3 x 10 repetitions, then left 3 x 10 repetitions.

**Note:**  
Exercise can be performed while sitting and standing.

**Accessory:**  
Resistance Band Forefoot Lifter (628174)

### S8 Foot Arch Stretcher 16002080

The arch of the foot is of great importance for a correct foot function since it acts like a shock absorber. The foot has a longitudinal arch and a cross arch. Thus, the body weight is mainly carried by three points: Heel, big toe joint and small toe joint. Foot problems such as flat-, skew- or splayfeet are often caused to a subsidence of the arch of the foot. However, a too pronounced arch also significantly impairs the function of the foot – we hereby speak of hollow foot.

Through this exercise the muscles are massaged and kneaded. This way the muscles keep smooth, soft and more flexible.

**Exercise scope:** At the beginning or depending on your personal feeling, stand on the half-spheres with more or less body weight and let the pressure point effect on your foot muscles. Initial pressure sensitivity is not uncommon when starting the practice; many people experience this, but it usually subsides after frequent use.

### S9 Plantar Stretcher 16002090

The plantar fascia (plantar aponeurosis) is a tendon plate found on the sole of the foot. It starts at the heel bone (calcaneus) and extends to the ends of the metatarsal bones (metatarsals). It is partly responsible for maintaining the upright longitudinal arch of the foot and prevents the bony arch from being pressed down under the weight of the body. If this tendon plate is chronically overloaded (=missing regeneration), it comes to an irritation and irradiating pain at the heel bone, so that foot problems are inevitable consequences. This is mostly expressed through stabbing pain when the foot is exposed to stain. Almost always various sports such as cycling, running or ballet—all those where much pressure is applied to the Plantar fascia or where the foot needs to be bent or stretched to a high degree—are the trigger. But also, workload e.g. long standing periods or incorrect footwear are main factors. The results are significant pressure pain often with swelling under the heel. The complaints in the heel can disappear after a few days without any specific therapy. But in case pain lasts for a longer period, it may result in chronic ailments. Persistent irritations cause micro-injuries of the tendon tissue that result in most cases in an ossification of the tendon.

A heel spur (calcaneal spur) is formed, which can be seen in the X-ray image as thorn-like bony protrusion from the bottom of the heel bone (calcaneal tuber). The therapy is very diverse. After reducing the burden, cold application is recommended (rubbing with ice cubes, cold foot baths). Furthermore, physical therapy exercises focusing on stretching the calf and foot muscles for mobilization are recommended. The exercise station Foot Workshop S9 Plantar Stretcher enables to stretch exactly this area also preemptively by exercises over: The first ray of the foot—Plantar fascia—Achilles tendon—calf muscles. This regenerative stretch is especially recommended after strenuous activity.

**Exercise scope:** Starting position (barefoot): Place the heel on the base plate and the foot in a relaxed position along the ramp. Then, move the knee forward to build up the stretch before increasing the intensity by carefully applying pressure to the foot. At this station, slow, careful practice and a balanced alternation of pressure and release are particularly important.

Dear customer,

the Pedalo products shown are made of wood, a natural product. Deviations in color, structure, grain, small knot inclusions or scars prove the genuineness of our product. They are natural and therefore a sign of the use of original wood and no reduction in quality.

#### Warnings:

Strictly observe! Use the device exclusively for the purposes described in these instructions (and any other Pedalo product information that may be possibly enclosed with the product and refer to this product).

Warning! Only use the device on level ground with even and firm surfaces. Protect item from damp and wet. After longer or frequent use, check that the screws are tight and retighten if necessary.

This device is meant for using from persons at the age of 14 up. Maximum weight capacity approx. 150 kg.

Copyright by Holz-Hoerz GmbH · Copies (even partial) not permitted.

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Discover more products for play, sport and  
therapy at [www.pedalo.de](http://www.pedalo.de)

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Copyright by Holz-Hoerz GmbH · Nachdruck auch auszugsweise nicht gestattet.