# pedalo®

## Ganzkörper-Trainingsprodukte Full-body training products

20050050, 18170038, 13029050, 18160038, 13006050,  
13007050, 18190038, 13008050, 18180038

[www.pedalo.de](http://www.pedalo.de)

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
This trainer was developed by us to strengthen the back, prevent falls, recover more quickly from injuries and alleviate existing complaints and pain. From children to seniors, everyone has fun training with it, because it is simple, efficient and enjoyable. The rotation disk made of natural birch multiplex serves to mobilize, relax and strengthen the deep back muscles and the lateral trunk muscles. Because it is precisely the rotational movement that can release vertebral blockages, tension and muscle pain from the sacrum up to the neck. In addition, rotational exercises around the entire body axis stimulate the sense of balance. The lateral abdominal and back muscles can be trained lying down. The knee rotators can be trained in the one-leg stand with lunge and the knee rotation is stabilized. Standing, the abductors are strengthened and the leg axis strengthened. The disc is highly stable and can be used on almost any floor without any problems.

**Specification:** Birch multiplex, clear lacquered. ø 38cm, standing height 5cm. With anti-slip feet. Weight capacity up to 100kg. 3kg. **18180038**

EN

Strengthen the back, prevent falls, reduce injuries, discomfort, and pain. Push motor development. The multifunctional Rotation Trainer combines twisting, circling (rotating) and through the accessories also rocking and circling on an inclined plane. The high-quality precision pivot bearing meets the requirements of physio practices, rehabilitation centers, sports performance centers, sports clubs, fitness centers and thus creates the conditions for successful training of body perception and body stability. The rotating movement relaxes, mobilizes, and strengthens the lateral torso and back muscles from the sacrum up to the neck. Fascial adhesions, vertebral blockages, tension, and muscle pain can be released. 360-degree rotation exercises around the complete body axis stimulate the sense of balance and stabilize the knee and hip rotators when standing on one leg. The disc has a high stability even on slippery floors. In combination with the accessories, the Rotation Trainer is a true all-rounder when it comes to sports, fitness, and health. Resistance stoppers/sliding resistors regulate the speed of rotation.

**Specification:** Birch multiplex, clear lacquered. ø 50cm, height 7cm. Grid holes at the top for the optional accessory Pedalo Footboards Seesaw (can be ordered under item no. 653101). With very high quality pivot bearing. 4 resistance stoppers with felt coating to regulate the rotation speed. 3 integrated magnetic holders on the underside for attaching the optional half-spheres in different positions (can be ordered under item no. 159072). Incl. rubberized anti-slip feet. Weight capacity up to 150kg. 6.7kg. **20050050**

### Accessory: Wood Half Sphere

The Rotation Trainer is equipped with magnetic holders for the wood half-spheres on the underside. By attaching two magnetic half-spheres next to each other (at a short distance apart), an inclined rotation disk is created, which has a dynamic effect during eccentric loading and accelerates the rotation around the body's longitudinal axis. If the two half-spheres are positioned opposite each other (on the central axis), even the slightest body imbalances will set the newly created Rotating-Balance Board into a seesawing and rotating motion.

**Specification:** Solid beechwood half-sphere, oiled, 6cm high, ø 10cm. Easy to attach thanks to the magnetic plug-in system. **159072**

### Accessory: Footboards Seesaw

Attachable foot rockers additionally increase the exercise movement requirements toward pronation and supination. The Footboards Seesaw can also be used individually as a training device to strengthen the ankle joints.

**Specification:** Birch multiplex, clear lacquered, 1 pair each 35x12x5.5cm. With fixing mechanism for attaching to Pedalo products that are equipped with suitable grid holes. Incl. small self-adhesive anti-slip strips for the underside and 4 larger anti-slip stickers for the top side. Weight capacity up to 120kg. 1.3kg. **653101**

Lying down flat → conscious activation of the rotation

Inclined plane → reactive training to the onset of rotation

Seesaw → reactive training on the onset of rotation plus establishing balance

The Balance Seesaw convinces through its functionality and very stable platform. In therapy and fitness training, the balance seesaw enables a variety of valuable stabilization exercises in standing, sitting and supported positions.

**Specification:** Birch multiplex, clear lacquered.

**Balance Seesaw 500:** ø 50cm, angle of inclination 8°, height 9cm. Grid holes for accessories. Incl. rubberized anti-slip feet. Weight capacity up to 150kg. 4.4kg. **13029050**

**Balance Seesaw 380:** ø 38cm, angle of inclination 12°. With anti-slip feet. Height 7cm. Weight capacity up to 100kg. 2.1kg. **18170038**

**Note on the exercises:** If exercises are marked as suitable for this product group, they can be performed with all variants.

left / right

forth / back

Trains the coordination, stabilizes the ankle and strengthens the muscles from the foot to the trunk. A Balance Top places greater demands on balance than a balance seesaw. The unique design with two boards allows it to be used on different floor surfaces while offering a high degree of stability.

**Specification:** Birch multiplex, clear lacquered.

**Sport+Kids:** ø 50cm, angle of inclination 13°, height 13cm. Grid holes for accessories. Incl. rubberized anti-slip feet. Weight capacity up to 120kg. 4.8kg. **13006050**

**Therapy+Fitness:** ø 50cm, angle of inclination 9°, height 9.5cm. Grid holes for accessories. Incl. rubberized anti-slip feet. Weight capacity up to 120kg. 4.6kg. **13007050**

**Balance Top 380:** ø 38cm, angle of inclination 12°. With anti-slip feet. Height 7cm. Weight capacity up to 100kg. 2kg. **18160038**

**Note on the exercises:** If exercises are marked as suitable for this product group, they can be performed with all variants.

The Trimm-Top is the balance, sport and therapy top made of wood for functional coordination training at the highest level. Ideal for one-legged balance exercises while standing, sitting and arm support.

**The benefits of use:**

- Effective balance and posture training
- Strengthening of foot, leg and back muscles
- Improvement of your reaction time
- Enables a range of motion up to the maximum possible angle of inclination
- Turned upside down, it is used to stretch and strengthen the cross arch of foot
- Highly functional spinner top, made of the natural material wood for a smooth movement pattern even under load
- For children, adults, seniors

**Specification:** Birch multiplex, clear lacquered.

**Trimm-Top 500 Vario:** ø 50cm, grid holes for accessories. Magnetic plug-in system on the underside for attaching the half-spheres as required (3 possible positions). 2 magnetic wooden half-spheres 6cm high made of beech. Stand height 7.5cm, angle of inclination 14°. Weight capacity up to 120kg. 2.1kg. **13008050**

**Trimm-Top 380:** ø 38cm, standing height approx. 7cm, angle of inclination 15°. Half-sphere made of solid beech, firmly mounted in the center of the wooden disc. Weight capacity up to 100kg. 1.15kg. **18190038**

**Note on the exercises:** If exercises are marked as suitable for this product group, they can be performed with all variants.

The possible positions of the half-spheres on the Rotation Trainer 500 and Trimm-Top 500 Vario are shown below (whereas with Balance Top or Seesaw, the positions/movement patterns are predefined through mounted rubber buffers):

**Exercise 1**

*Possible with: Rotation Trainer ▪ Twister*

Twist side to side (without the resistance stoppers). Your upper body rotates against your hips.

**Variation for the Rotation Trainer:** Twist against the resistance of 1 to 4 stoppers.

**Exercise 2**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

**Starting position:** Stand lengthwise to the half-spheres.

Try to stabilize the training device from a stable, two-footed stance.

Rock side to side, maintaining your directional position.



**Exercise 3**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

**Starting position:** Stand crosswise to the half-spheres.

Try to stabilize the training device from a stable, two-footed stance.

Rock back and forth.

**Exercise 4**




*Possible with: Rotation Trainer*

**Starting position:** Stand lengthwise to the half-spheres.

Try to stabilize the training device from a stable, two-footed stance.

Perform a controlled ¼, ½, ¾, or full rotation and stabilize in the starting position.

Do this exercise in both directions.

**Exercise 5**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top ▪ Twister*

Perform without using the stoppers and half-spheres (Rotation Trainer).

For this exercise, anchor the band to the wall bars or have a partner hold it.

Pull the resistance band with both hands parallel to your body, at chest height to the body, over your head, or diagonally past your body with both hands.

**Variation for the Rotation Trainer:** Repeat with the half-spheres (no stoppers).

**Exercise 6**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top ▪ Twister*

Perform without using the stoppers and half-spheres (Rotation Trainer).

Pull the resistance band, coming from behind, with both hands parallel to/past your body, at chest height past the body, over your head (throw-in), or diagonally past your body with both hands. Maintain a steady posture – do not rotate. Counteract as needed.

**Variation for the Rotation Trainer:** Repeat with the half-spheres (no stoppers).

**Exercise 7**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

Lean your upper body forward and shift your weight onto your arms.

As before, additionally rock from side to side.

Repeat both parts of the exercise, this time rocking back and forth. To do this, first turn the training device 90 degrees so that the half-spheres are aligned according to your desired direction of movement.

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

**Exercise 8 (as in 7 but with knees raised)**

*Possible with: Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

Lean your upper body forward and shift your weight onto your arms. Hover your knees just above the floor.

As before, additionally rock from side to side.

Repeat both parts of the exercise, this time rocking back and forth. (Half-spheres/training device aligned according to the other desired direction of movement.)

**Variation for the Rotation Trainer:** Perform without using the stoppers and repeat the exercise.

### Exercise 9

**Possible with:** *Rotation Trainer* ▪ *Trimm-Top* ▪ *Twister*

**For the Rotation Trainer:** Remove the stoppers and place the half-spheres front-to-back (relative to your position) to allow for side-to-side movement.  
**For the Trimm-Top Vario:** Place 1 half-sphere in the center.

Lean your upper body forward and shift your weight onto your arms. Hover your knees just above the floor.

Rotate the coordination device to the right and left. Try to keep it level while doing so.

**Variation for the Rotation Trainer:** In the starting position, the half-spheres are positioned crosswise to you.

### Exercise 10

**Possible with:** *Rotation Trainer* ▪ *Balance Seesaw* ▪ *Balance Top* ▪ *Trimm-Top*

Lean your upper body forward and shift your weight onto your arms. Hover your knees just above the floor.

Try to keep the training device level and extend one leg.

Repeat the exercise, this time with the training device in a changed position (half-spheres now crosswise/device turned by 90 degrees).

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 11

**Possible with:** *Rotation Trainer* ▪ *Balance Seesaw* ▪ *Balance Top* ▪ *Trimm-Top*

Stabilize the coordination device with one hand. Then extend the other hand and, diagonally, the leg.

Repeat the exercise using a different device placement (half-spheres front-to-back relative to you/side-to-side movement).

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 12

**Possible with:** *Rotation Trainer* ▪ *Balance Seesaw* ▪ *Balance Top* ▪ *Trimm-Top*

Get into a push-up position.

Hold the coordination device level while alternating leg extensions.

As before, but this time alternate pulling each knee toward your chest.

Repeat the exercise using a different device placement (half-spheres crosswise/device turned by 90 degrees).

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 13

**Possible with:** *Rotation Trainer* ▪ *Balance Seesaw* ▪ *Balance Top* ▪ *Trimm-Top*

Get into a push-up position and lower your body.

Keep the training device level while extending one leg.

Repeat the first exercise using a different device placement (half-spheres crosswise/device turned by 90 degrees).

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 14

**Possible with:** *Rotation Trainer* ▪ *Trimm-Top* ▪ *Twister*

**For the Rotation Trainer:** Remove the stoppers and place the half-spheres front-to-back (relative to your position) to allow for side-to-side movement.  
**For the Trimm-Top Vario:** Place 1 half-sphere in the center.

Get into a deep push-up position and keep the coordination device level.

Then rotate the training device to the right and left.

**Variation for the Rotation Trainer:** In the starting position, the half-spheres are positioned crosswise to you.

### Exercise 15

**Possible with:** *Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top ▪ Twister*

Get into a forearm support position and stabilize the coordination device with your feet on it.

While doing so, keep your pelvis and the training device as level as possible.

As before, but now alternately extend one hand forward.

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 16

**Possible with:** *Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

Get into a push-up position and stabilize the coordination device with your feet on it.

Lower yourself into a deep push-up position and slowly rock the board from side to side several times with your feet.

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

### Exercise 17

**Possible with:** *Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top ▪ Twister*

From a supine position, push your hips upward, lifting your hips as high as possible.

Hold this position, maintain your pelvis pushed upwards and the training device level.

As before, but place only one foot in the center and extend the other leg upward.

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise. Additionally rotate the twist board with your feet to the left and right.

### Exercise 18

**Possible with:** *Rotation Trainer ▪ Twister*

**For the Rotation Trainer:** Remove the stoppers.

From a reverse plank, walk your feet in a full circle around the coordination device.

Keep your hips high and the board level as you perform the movement.

### Exercise 19

**Possible with:** *Rotation Trainer ▪ Balance Seesaw ▪ Balance Top ▪ Trimm-Top*

From a reverse plank, alternately lift one leg while stabilizing the board's side-to-side motion. Keep your hips high and as steady as possible throughout.

As before, but with the half-spheres positioned to allow for front-to-back rocking.

**Variation for the Rotation Trainer:** Remove the stoppers and repeat the exercise.

Dear customer,

the Pedalo products shown are made of wood, a natural product. Deviations in color, structure, grain, small knot inclusions or scars prove the genuineness of our product. They are natural and therefore a sign of the use of original wood and no reduction in quality.

#### Warnings:

Warning! The device is to be used with caution since it requires great skill to avoid falls or collisions causing injury to the user or third parties. In case of doubt first use the device only with assistance and/or appropriate protective equipment. Use only under adult supervision.

Only use the device on level ground with even and firm surfaces. Protect item from damp and wet. After longer or frequent use, check that the screws are tight and retighten if necessary.

**Additional note on the Rotation Trainer:** Caution! When using with the optional Wood Half-Spheres/Footboards Seesaw, the Rotation Trainer becomes a very challenging exercise equipment.

A particular degree of caution, attention, circumspection, and a high level of motor skills are required, especially during initial attempts and/or with standing exercises.

The devices depend on the user's coordination skills. The smaller 380 devices, featuring a double-deck design, are recommended for ages 3 and up, while the Trimm-Top 380 is recommended for ages 5 and up. However, the larger devices in the 500 series, Balance Seesaw, Balance Top Sport+Kids, and Trimm-Top Vario, are recommended for children aged 7 and above, while Rotation Trainer and Balance Top Therapy+Fitness are recommended for children aged 14 and above. The maximum weight capacity ranges from 100 to 150 kg, depending on the model.

Copyright by Holz-Hoerz GmbH · Copies (even partial) not permitted.

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Discover more products for play, sport and  
therapy at [www.pedalo.de](http://www.pedalo.de)

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Copyright by Holz-Hoerz GmbH · Nachdruck auch auszugsweise nicht gestattet.