pedalo®


Fußwerkstatt
Foot Workshop

16002000

www.pedalo.de

Ossa digiti pedis Zehenknochen Phalanges	Phalanges proximales Zehengrundglieder Proximal phalanges
Phalanges distales Zehenendglieder Distal phalanges	Metatarsalia Mittelfußknochen Metatarsal bone
Phalanges mediae Zehenmittelglieder Middle phalanges	Ossa tarsi Fußwurzelknochen Tarsal bone

1. Talus Sprungbein Ankle bone	4. Os cuneiforme I Keilbein Sphenoid bone	
2. Calcaneus Fersenbein Heel bone	5. Os cuneiforme II Keilbein Sphenoid bone	
3. Os naviculare Kahnbein Scaphoid bone	6. Os cuneiforme III Keilbein Sphenoid bone	7. Os cuboideum Würfelbein Cuboid bone

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
Our feet are the fundamental base of our body whose importance we underestimate. The human foot is a fantastic construction of 26 bones, from many ligaments, tendons and nerves. If healthy it carries our full body weight painlessly and perfectly cushions each step gently. It regulates the balance and provides a wide range of movement possibilities. Only with healthy feet, we have the necessary and important security of movement from a young age, which we need throughout our life for fall prevention - reactively but also preventively. Only with healthy feet, we have fun in the movement and are able

to stand, walk, run, jump, dance and bounce. To achieve all these tasks, the foot relies on a healthy and resilient vault structure. A majority of the population is suffering from foot problems of which two-thirds are often caused by (wrong) movement habits and inappropriate biomechanical stress on muscles and joints, what leads to skew-, flat- or splayfeet. This inevitably means an over-stressing of foot, knee and hip joints and spine. Often this is the beginning of a chain of larger foot problems and in its further course it can lead to chronic inflammation and increased signs of wear and tear

in the joints, associated with pain. Shoe inlays, prescribed by a physician, support the foot position and may relieve pain but they also place the foot in its function passively and do not contribute to strengthen the insufficient foot musculature. We from pedalo therefore recommend you a daily training dosage for your feet. This is the only way to strengthen the arch of the foot and even to improve incipient foot deformities. Foot training begins with bare-foot walking. Then it includes a daily mobilization of feet, associated with an active muscular structure workout.

In everyday life and in sports healthy feet are a blessing and the basis for pain-free movement. Almost all malformations of feet are due to muscular deficiencies or imbalance of the foot and leg muscles – and this often already from an early age on. To strengthen the muscles a continuous training is essential. The nine exercise stations of the Pedalo-Foot Workshop are convincing through their versatility and activate all systems and structures, which are important for a healthy and functional foot. Thus the different stations are intended to activate the feet, sensitize and massage the foot soles, stabilize the ankle joints, strengthen the foot and lower leg musculature as well as to stretch the foot arch and the plantar fascia. The Pedalo-Foot Workshop is suited ideally as training circuit for prevention and rehabilitation. It is used in fitness & health studios, podiatry, physiotherapy and rehabilitation facilities, in retirement homes or even in wellness hotels but also in the course of operational health management. The shelving frame offers a space-saving and clear possibility for storage.



- S9 *Plantar stretcher*
- S1 *Twister*
- S2 *Toe grip*
- S3 *Senso*
- S4 *Seesaw*
- S5 *Top*
- S7 *Forefoot lifter (New version)*
- S6 *Heel lifter*
- S8 *Foot arch stretcher*

S1 Twister 16002010

When doing the twist we do a dynamic movement and almost all muscles of the foot and the lower leg are involved. As a starting position place both feet side by side. Now the heel is lifted, turned on the forefoot to side and lowered down again. Then the forefoot is lifted, turned on the heel side and put down again. Repeat the same exercise in opposite direction now. Perform the exercise in a slow and controlled manner.

Exercise scope: Repeat 10 times left and right.

S2 Toe grip 16002020

The toe grip strengthens and relaxes the foot and the toe muscles. Similar to a caterpillar, pull the foot forwards by stretching and tightening of the toes. The cross grooves support the forward movement. Make sure that all toes grab in the grooves. Then carry out the contrary toe movement. To this extend the toes wide upwards.

Exercise scope: Repeat 5 times.

S3 Senso 16002030

The Senso board is like an acupuncture on the sole of the foot and takes positive influence on the nervous system and the activation of the metabolism. The foot muscles are stimulated and relaxed at the same time.

Exercise scope: Stand on the wooden nubs for 3 minutes but change the foot position several times.

Exercises with the Foot Workshop

S4 Seesaw 16002040

Footboards Seesaw offer a linear movement of the foot into pronation and supination (strain of inner/outer side of the foot). Beside of the foot muscles, the stabilizing musculature of the lower leg is activated. The combination of these muscles affects our posture stability.

Exercise scope: Repeat 10 times parallel movement left-right, then stabilize for 10 seconds in horizontal position.

S5 Top 16002050

Footboards Top tilt on all sides and are very close to the function of the ankle joint and its mobility. The requirements to the coordination skills are here higher than with the Footboards Seesaw. Therefore this station serves for mobilization and the stabilization of the ankle joint.

Exercise scope: Train the ankle joint by doing several alternating circular movements left and right (clockwise, counterclockwise or both at the same time) against the resistance of the Footboards Top. Stabilize the Footboards Top out of motion and freeze the position for 10 seconds. Stand on one leg right and left alternatingly.

S6 Heel lifter 16002060

Exercises on the Heel lifter activate the interaction of foot muscles with the lower leg muscles and give the ankle joint its necessary stability. Already the starting position - forefoot is placed on the round wood, heel on the grey plate - means a stretching in the area of Achilles tendon/calf muscles. By lifting the heel in the horizontal position the muscles around the foot, ankle and lower leg, as well as the sense of balance and responsiveness are particularly required. The lifting and lowering of the heel should be executed slowly and evenly.

Exercise scope: Repeat 10 times.

Exercises with the Foot Workshop

S7 Forefoot lifter 16002070

People with foot weakness must have much concentration not to fail in taking steps. Every little bump means a danger zone. So now, the normally almost automated movement process of walking, represents for the person concerned - in the truest sense of the word - a task, that he must perform step by step attentively. Unfortunately this always affects his perception of the environment adversely. Very often people suffering from weak foot dorsiflexion, place their foot not with the heel first but with the forefoot or the entire sole of the foot. A natural process of rolling of the foot and forward swing of the leg is no longer possible for foot weak dorsiflexion. The person concerned becomes cramped, walks abnormally and sometimes also even with clenched toes. Swinging the leg forward is done by a circular hip movement (Circumduction). In the long term people moving themselves such way, complain about pain in hip and pelvis. Anatomy: The muscle tibialis anterior supports the dorsal flexion, respectively extension, supination, and adduction of the foot. In addition, his tension stabilizes the ankle joint, especially in the first phase of ground contact when walking or running. He is the antagonist of the muscle 'peroneus longus'. The paralysis of the muscle 'tibialis anterior' leads to the so called steppage gait. With the Forefoot lifter you activate the muscles responsible for a better gait through targeted exercises.

Exercise: Starting position: Training without shoes. You can choose working against 1-4 resistance bands (the more bands the more resistance); the forefoot is placed that way so that the resistance band runs diagonally across the basic joints of the toes (toes look beyond the resistance band); heel is placed on the half round wood. Intensity: Depending on the degree of movement, the knee of the exercise foot is positioned vertically above or behind the heel. Note: >> Maximum degree of movement or movement amplitude = forefoot is stretched to its maximum and the knee is behind the heel as far as possible. >> Maximum tension = Starting position knee over heel and the forefoot needs to raise to its maximum. The leg itself remains fully relaxed during the exercise execution. The concentration is the front of the foot and anterior tibial muscle. The toes should remain fully stretched during exercise and are not pulled up.

Exercise scope: 3 times 10 repetitions per foot. It is important to keep up the number of repetitions, thus please select the number of resistance bands accordingly so that the exercise series can be performed well. Do not

change exercise foot during the individual exercise series. It means: At first with the right foot: 3 x 10 repetitions. Only then with the left foot 3 x 10 repetitions.

Notice:
Exercise can be performed while sitting and standing.

Accessory:
Resistance Strap Forefoot lifter (628174).

S8 Foot arch stretcher 16002080

The arch of the foot is of great importance for a correct foot function since it acts like a shock absorber. Foot problems such as flat-, skew- or splayfeet are often caused to a subsidence of the arch of the foot. However, a too pronounced arch also significantly impairs the function of the foot – we hereby speak of hollow foot. The foot has a longitudinal arch and a cross arch. Thus the body weight is mainly carried by three points: Heel, big toe joint and small toe joint. Depending on the strain of the foot, the Pedalo Foot arch stretcher exerts a stronger, punctual stimulus on the muscles. Through this exercise the muscles are massaged and kneaded. This way the muscles keep smooth, soft and more flexible. At the beginning or depending on your personal feeling, stand on the hemispheres with more or less body weight and let the pressure point effect on your foot muscles. Initial pressure sensitivity is not uncommon at the start of the practice, many people do have so, but it is usually dissolved after frequent use.

Exercise scope: Let a pressure point act on your sole of the foot. Practice alternately or with both feet together for 10 seconds. Repeat 3 times.

S9 Plantar stretcher 16002090

The tendon plate Plantar fascia (plantar aponeurosis) is found on the sole of the foot, where it starts at the heel bone (tuber calcanei) and continues to the end of the metatarsal bones (Ossa metatarsals). It is jointly responsible for the maintenance of an upright position of the longitudinal arch of the foot and prevents, the bony foot arch to be pressed the way down under the body weight. If this tendon plate is chronically overloaded (=missing regeneration), it comes to an irritation and irradiating pain at the heel bone, so that foot problems are inevitable consequences. This is mostly expressed through stabbing pain when stain is placed on the foot. Almost always various sports such as cycling, running or ballet – all those where much pressure is applied to the Plantar fascia or where the foot needs to be bent or stretched to a high degree – are the trigger. But also work load e.g. long standing periods or incorrect footwear are main factors. The results are significant pressure pain also with swelling under the heel. The complaints in the heel can disappear after a few days without any specific therapy. But in case pain lasts for a longer period, it may result in chronic ailments. Persistent irritations cause micro-injuries of the tendon tissue that result in most cases in an ossification of the tendon. A heel

spur (calcaneal spur) is formed, which can be seen in the X-ray image as thorn-like bony protrusion from the bottom of the heel bone (calcaneal tuber). The therapy is very diverse. After reducing the burden, cold application is recommended (rubbing with ice cubes, cold foot baths). Furthermore physical therapy exercises with stretching the calf and foot muscles for mobilization are recommended. The exercise station 'Foot Workshop S9 Plantar stretcher', enables the practitioner to stretch exactly this area also preemptively by exercises over: The first ray of the foot - Plantar fascia - Achilles tendon - calf muscles. Especially after heavy strain it is recommended to do this regenerative stretching exercise.

Exercise scope: Starting position (barefoot): The heel is on the base plate and the foot is in relaxed position along the ramp. Execution: The knee is moved forward so as to obtain the muscle elongation. The exercise is carried out under load of the foot whereby the pressure is increased slowly. Practice as follows: Stretch slowly – soften the pressure – stretch little harder – soften the pressure – and so on.

Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore signs of using original wood and no quality reduction.

Please note:

WARNING! The Pedalo-Foot Workshop is only to be used for the purposes described in this user manual.

WARNING! Place the unit only on flat level ground.

After longer or frequent use, check that the screws and all installed parts are tight and retighten if necessary. Defective and worn parts should be replaced by original replacement parts immediately.

Protect item for damp and wet. This device is meant for using from persons at the age of 6 up. Maximum weight capacity approx. 150 kg

Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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