pedalo®

Pedasan 220/320

18250022, 18250032

www.pedalo.de

18250022

18250032

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EN

Pedasan - a valuable sport and training device for healthy feet and intact balance.

Pedasan in sports and therapy

Exercising barefoot on Pedasan helps to strengthen the foot musculature and to boost the neglected foot vault. The blood circulation in feet and legs will be stimulated. With shoes the Pedasan can be used also for coordination exercises which support the ability of balance and reaction as well as an upright body posture. At the big Pedasan 320 the pressure of the feet spread evenly to a bigger surface. Therefore, the feet have more grip and the exercises are easier. Through the higher stand, users need more courage as at the smaller Pedasan 220.

Pedasan in play

Movements on the Pedasan have a very high demanding character, contribute to the enrichment of the sports offer in schools and clubs, and are also a valuable aid for testing and perfecting the ability to balance. After previous balancing experiences on benches, balance beams, etc., the use is recommended in physical education classes, break time sports, recreational activities, club groups, occupational therapy, and psychomotor movement promotion.


Specification: Side parts made of birch multiplex, with rubber tires. Beech wood round bars. Width 38cm. Max. weight approx. 100kg.

18250022 ø 22cm. 2.3kg

Unit

18250032 ø 32cm. 3.9kg

Unit

Accessory: Support Rope (556001-01)

Concentrated and relaxed walking and standing on the Pedasan is enough to train the foot muscles.

Also usable as a stool.

All exercises are practicable with the Pedasan 220 and 320.

The training improves the vestibulomotor skills and the balance. If you train without shoes, your foot arch is built effectively. Pedasan 320 has a lower pressure on the foot soles, because it has a larger radius than the 220.

Exercise while standing

Stand with both feet on the Pedasan. Keep the balance. Slowly roll forward and backward. Helpful at the beginning: Walk on the Pedasan along a wall.

Variation: Stand with one feet on the Pedasan. The other one stays on the ground. Roll the Pedasan forth and back.

Exercise while walking

Stand with both feet on the Pedasan. Slowly walk forward and backward. Helpful at the beginning: walk along a wall on the Pedasan.

Variations:

- Change viewing direction (forward, to the left, to the right, looking up, close your eyes)
- Move on it in a squatting position
- Swing your arms

Exercises for self-massage and regeneration

The Pedasan is also a perfect training product for self-massage and for relaxation of the connective tissue. Roll with the Pedasan on different muscle groups. At painful or tensed points, hold the position for some seconds and then continue. By supporting yourself with your arms and inactive leg, you vary the pressure on the Pedasan and at the same time increase the stability of your body.

Exercises in the support

Roll the Pedasan forward away from the body and backward to the body.

Variation: Roll forward and backward with your feet on the Pedasan.

Exercises with two Pedasan

Put your feet and hands each on one Pedasan. Then:

- roll forward/backward
- roll the Pedasan towards/away from the other

Variations:

- Each foot is on one Pedasan. Now roll forward/backward on different heights.
- Stand on a Pedasan and start moving. Change to the other Pedasan and continue rolling.

Exercises with the ball

Always go forward/backward in rhythm. In doing so:

- bounce a ball (left/right hand)
- make a throw-in over your head against a wall
- throw the ball into a basketball hoop, goal, box, hoop...
- circle the ball around your hips or pass it through your legs

Exercises with ball and partner

Always go forward/backward in rhythm. In doing so:

- throw/catch a ball directly to/from the partner
- throw/catch a ball indirectly through the ground
- throw/catch with one hand
- make a throw-in over your head
- head the ball back
- try throwing techniques from volleyball, such as digging

Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of the material. They are natural and therefore features of using original wood and no quality reduction.

Please note:

The Pedalo Pedasan is only to be used for the purposes described in this user manual.

Warning! The device is to be used with caution, since it requires great skill, so as to avoid falls or collisions causing injury to the user or third parties. In case of doubt first use the device only with assistance and/or appropriate protective equipment. When used by children the supervision by an adult is recommended.

Caution! Place the unit only on flat level ground. Protect item for damp and wet.

After longer or frequent use, check that the screws are tight and retighten if necessary.

This device is meant for using from persons at the age of 5 up.

Maximum weight capacity approx. 100 kg.

Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

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