pedalo®

Fußmassage-Regenerationsmatte Foot Massage Regeneration Mat

143028, 23143028

www.pedalo.de

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
Promotes blood circulation when standing for long time.

Water resistant and adheres to the shower tray. Regenerate during the shower.

Prevention through sensory stimulation before sports. Fast regeneration after sports.

- promotes blood circulation with venous disorders and cold feet
- Tension and relaxation after and during long sitting and standing at work
- improved regeneration after sports activities
- Revitalization for tired and heavy legs
- Reflex stimulation of the whole body
- Stimulation of the whole body in case of sensitivity disorders
- waterproof: ideal in the shower at home and in sports


Our feet are the fundamental base of our body and their significance should not be underestimated.

It is generally known that our feet reflex our whole body. Our feet have a very complex structure with many bones, tendons and muscles that are created in a wonderful way.

Feet carry us for a lifetime and are still too often neglected. In its fine and highly complex construction, they are unique and do not require any improvement. In their construction they are not equal with anything. Unfortunately, we often notice them only when they call for help. From the very beginning of our ability

to stand and walk, the entire foot structure strengthens so that it is constantly available to us. Walking barefoot is now mostly no longer possible anywhere. But it can be replaced by standing or walking on pebbles. In this way it can be seen how sensitive our feet are. Regular daily exercises of this type comply with reflexology massage. It is a health care of high quality and also ensures the health of our entire body. It depends on the mindful treading, consciously and with feeling, so that the feet become perceptible again.

Recommended by
Brigitte Wilhelmer,
foot reflexology practice

General Notes: This product is made of natural and synthetic rubber. Slight discoloration of the material is natural in the production of products made from natural rubber and therefore no reason for complaint.

Sensory system and senior age

As a result of today's lifestyle on hard and straight ground conditions, differentiated and especially natural surface stimuli are lacking. The sensory systems of the soles of the feet atrophy, changing balance information and corresponding neuromuscular reactions and adaptations of the total muscle system (dysfunctions of the foot and the ascending muscle chains) are missing, what may lead to dysfunctions of the foot and the ascending muscle chains.

Special added value for seniors or people who depend on being visited by caregivers: Many get their legs wrapped, bandaged or wear compression stockings. Unfortunately, this can result in the foot being deprived of its natural task of circulating blood and decongesting veins. Therefore, the hint: before wrapping/bandaging and while wearing compression stockings, the foot massage regeneration mat (for example, during morning hygiene) provides positive and effective stimulation of the body's natural functions.

A foot sole massage is not only a benefit for the feet, but above all has a significant effect on health and well-being. Due to the spherical curvature of the Pedalo Foot Massage Regeneration Mat, all stimulus points of the soles of the feet are stimulated simultaneously. Just 3 minutes a day is enough not only to release your foot musculature and consequently the deep musculature of the whole body, but also to strengthen the cross arch of foot and to prevent foot suffering. Blood circulation and purification of the legs are stimulated and activated. It gives a pleasant feeling of warmth. In sports it is used for quick recovery from stress. The Massage Mat is ideal to quickly take and use anywhere, for example, when traveling, after a workout or in the shower.

Specification: Made of natural and synthetic rubber, ø 25 cm. Maximum load up to 150 kg. 0.6 kg.

143028 Blue

23143028 Anthracite

Notes for use

The talc on the top protects the mat. It can be practiced on it without hesitation. If it bothers you visually, it can be washed off with washing-up liquid and water. The foot regeneration and massage mat does not replace walking barefoot.

How to clean the Massage Mat

Simply hold under water and brush off. You can also put it in the dishwasher at a low temperature setting.

1. Start with the right foot. Place it on the massage mat so that the outside of the foot is towards the center of the mat.
2. Then switch to the left foot and repeat the exercise.
3. Place both feet parallel on the massage mat. The inner sides of the feet should enclose the center of the massage mat.
4. Sit on a chair and place your feet on the mat. Apply alternating pressure on the mat with the feet. Vary the pressure on the soles of the feet.
5. Extend the exercise by alternately lifting your feet off and then letting them dip back into the nubs.
6. Stand upright and repeat exercises 4 and 5.

Notes:

- Do each exercise for 1 minute at a time. Make sure your knees are slightly bent, toes are loose, and ankles are relaxed.
- It is normal to feel initial pressure pain. Stay relaxed and let the pimples effect on the sole of the foot. If you feel discomfort, it is advisable to extend the pauses or vary the exercise (shorten or return to a sitting/walking position instead of standing) or stop for the moment and rest. The body must first adapt to many stimuli step by step.
- In sitting position the feet can be placed on the mat as long as it feels comfortable.
- If you wish, you can also sit on the massage mat. This can help relax your neck, shoulders, arms and hands.

Breaks: We recommend taking a break of at least 30-60 seconds between exercises. If necessary, the pauses can be extended.

Bear in mind: If you feel pain the first few times using the mat, it's a sign that your feet need help. The causes are muscle and fascia tension, as well as deposits – it is not because of the mat. Persevere; after several uses, you will realize how good it is for you.

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Discover more products for play, sport and
therapy at www.pedalo.de

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Holz-Hoerz GmbH