pedalo®

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Our feet are the fundamental base of our body whose importance we mostly underestimate. The human foot is a fantastic construction of 26 bones, from many ligaments, tendons and nerves. The physician Dr. Christian Larsen discovered that nature has built our feet according to the spiral principle. The back foot (heel) rotates outward, the fore foot turns oppositely inward. Therefore, the so important lengthwise and crosswise arch are arised.

A majority of the population is suffering from foot problems of which two-thirds are often caused by (wrong) movement habits and inappropriate biomechanical stress on muscles and joints, what leads to skew-, flat- or splayfeet. This inevitably means an over-stressing of foot, knee and hip joints and spine. Often this is the beginning of a chain of larger foot problems and in its further course it can lead to chronic inflammation and increased signs of wear

and tear in the joints, associated with pain. Therefore, the feet get a great significance for the whole locomotor system as basis of the straightening and body posture. In its totality, the sole of the foot forms a sensory contact zone between body and earth. By the important perceptual signals about the nature and position of the subsoil, which it sends to the brain, the balance of the inner ear is supplemented. Also the horizontal position of the head is regulated and the sensomotoric balance of the body adjusted.

Optimal longitudinal and transverse arches of the foot are essential for a healthy foot. A healthy longitudinal and transverse arch construction provides a unique and inherently resilient construction, which allows for a statodynamic pressure distribution.

In addition to the numerous foot muscles, the structures of the plantar fascia and the long soles of the

foot are important factors in the development of numerous short muscles. In addition, the interaction of the structures of the upper and lower ankle joints allows the deflection of vertical and horizontal forces acting on the foot. Further, the ankle joints ensure the mobility of the foot skeleton against the lower leg. Decisive for the function and stability of the hindfoot is the active joint stabilization and control of the heel contact during the appearance by the musculature. The round basic shape and relative small area of the heel requires the immediate balanced activation of the stabilizing lower leg and foot muscles.

With the following Pedalo-products you can train these structures and characteristics which are decisive for a healthy foot.

Foot Torsion Trainer (14008000)

Our variable artist for your feet.

The unique foot torsion trainer perfectly adapts to your foot position. Through the individual segments you can see if your foot has a correct position or you train through them selectively the individual areas, which are considered together as a foot screw. Your foot should show an even spiral position from forefoot to hindfoot. Weak points, malpositions and thus instabilities become visible for you and can be trained in a unique way. With the help of two integrated fixation rods you get the strength, coordination and control from segment to segment decoupled / coupled with each other and can thus train even more specifically.

Specifications: solid beech. 35x12x5 cm. Maximum load up to 120 kg. 2 fixing rods. 0.2 kg.

Exercises with the Foot Torsion Trainer

Standing stabil on one foot

Now we practice the one-leg stand.

1. Lift one foot and remain stable. You should not move your upper body too much. Hold this for 20 seconds. Pause for 20 seconds. During this time, you can practice the other foot.
2. Now lift the leg that is in the desire and go into a standing balance. Hold the position as long as you can. Repeat this 3 times.
3. Describe a circle around yourself with the floating leg. The radius should be as large as possible. This is your goal. You can repeat this as many times as you enjoy and can perform it. Don't forget to practice clockwise and counterclockwise. Change your foot during the break.

Stable forefoot step position

1. Tand with the forefoot (reference bunion of the big toe and metatarsal bones are at the highest point) on the torsion trainer. Take the same position with the other foot on the other trainer. The soles of the feet are parallel to the floor.
2. Start weight shifting forward/backward, left/right, circling. Close your eyes while doing this, if you are already more confident.
3. Stand still and shift your weight to the front foot. Go heel to toe (no contact) and back up (as high as you can). Foot change. Repeat 20 times and give yourself a break of 30 seconds. Series: 4x.
4. Rotate your upper body to the left and right and perform the previous exercises.

Exercises with the Foot Torsion Trainer

Stabile basic position

1. Stand with the heel bone (calcaneus) on the felt highest point.
2. Walk on the spot and try to return to the starting position again and again.
3. Now perform a squat. Make sure that the knees do not bend inward or outward. Repeat the squat until the control over the movement decreases and the movements become increasingly impure. - Give yourself a break of 60 seconds.
4. Repeat the exercise 4 times.

Stable stance leg phase - single leg

1. Lift the back foot and perform the previous exercises again. Do not forget to change legs.
2. Go with the heel on the floor (with contact) and shift your weight forward. This way you can stretch and relieve the plantar side of the foot.

TIP: Vary the order between exercises and try never to touch the floor. The pauses are made on the product.

You can also close your eyes during the exercises.

Stable forefoot - double legged

1. Try to maintain the position by keeping the soles of your feet parallel to the floor.
2. Go with the heel to the floor and up again. Raise the heel as high as possible. Repeat this as many times as possible. Pause for 60 seconds and repeat 4x.
3. From the basic position of 1. Go into a squat and hold the position. Vary in the highs and lows. Important: a) hold the position b) dynamically go up and down. Repetitions 20x. Rest 60 seconds, series 4x.

Unstable basic position fixed

Unstable basic position one-legged

Unstable basic position free moving

Unstable basic position one-legged free moving

Unstable basic position

(free moving=remove wooden stick)

1. Familiarize yourself with the basic position and remain standing. Move your upper body to the left and right, rotate to the left and right and make sure your knees remain stable.
2. Alternately -slowly- lift the left and right foot.
3. Now combine 1. and 2.
4. Perform all movements slowly. Take 3 minutes for the complete sequence. Take a break of 60 seconds. Do it again, only 1x.


Tip: For leg stability we use the following: Take a resistance band and attach it to the door and go into the loop with one leg or hold it. Here you can let off steam and do strengthening exercises. Your advantage: The Foot Torsion Trainer in this combination will clearly show how wide your feet are and train them at the same time.

Unstable forefoot & Cross stance

(free moving=remove wooden stick)

1. Go onto the trainer in a controlled manner and hold the position for 30 seconds. As an increase you can close your eyes (no time limit).
2. Take out the bar and repeat exercise 1.
3. Alternately lift one foot slowly. You can pull the knee up to the nose.
4. Step onto a trainer with both feet and do a (target: deep) squat. Repeat as many times as you can. Rest 30 seconds. Repeat 4 times.
5. Start walking back and forth. Then turn right around so that you are going from bar to bar (trainer to trainer) in a circle and your body is turning completely 360° in a circle.

Tip: For leg stability, we use the following: Take a resistance band and attach it to the door and go into the loop with one leg or hold it. Here you can let off steam and do strengthening exercises. Your advantage: The Foot Torsion Trainer in this combination will clearly show how wide your feet are and train them at the same time.

fixed on both legs


free movement on both legs


one-legged fixed

one-legged free movable

Cross stance fixed


Cross stance free movable

Stable forefoot step position & Spiraldynamic

1. Stand in rope dancing position (in a line). Hold the position and begin to turn yourself, squat or close your eyes. Perform the movement in a controlled manner and very slowly.
2. Begin by shifting weight forward/backward and then move into a standing balance. Hold the position for 30 seconds. If you are already more confident, close your eyes.
3. Turn the front torsion trainer over and position yourself as follows:

Stable forefoot step position

Exercises with the Foot Torsion Trainer

4. Hold the position and lift the other leg into the air. Now begin to actively move the forefoot by pressing down on the outer edge of the forefoot. Then the ball of the big toe moves towards the floor.
5. Perform exercise 4. and begin to describe a circle around you. The radius should be as large as possible. This is your goal. You can repeat this as many times as you enjoy and can perform it. Don't forget to practice clockwise and counterclockwise. Change your foot when you pause.

Tip: For leg stability we use the following: Take a resistance band and attach it to the door and walk with one leg in the loop or hold it. Here you can let off steam and do strengthening exercises. Your advantage: the Foot Torsion Trainer in this combination will clearly show how wide your feet are and train them at the same time.

Spiraldynamic heel fixed

Dear customer,

the Pedalo-product you have chosen is made of natural wood. Differences in color, structure, grain, small included knots or scars, prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Please note:

Warning! The Pedalo-Foot Trainer is only to be used for the purposes described in this user manual. The device is to be used with caution, since it requires great skill, so as to avoid falls or collisions causing injury to the user or third parties. In case of doubt first use the device only with assistance and/or appropriate protective equipment. When used by children the supervision by an adult is recommended.

Warning! Place the unit only on flat level ground. Can be wiped off with a damp cloth. Do not use any solvent-based, acidic or abrasive detergents. Protect item for damp and wet.

This device is meant for using from persons at the age of 3 up.

Maximum weight capacity approx. 150 kg.

Copyright by Holz-Hoerz GmbH,
printed copies also partially are not allowed.

pedalo[®]

Academy – since 1963

EN

Your competence academy for
prevention, rehabilitation and increase performance.
Concepts, coaching and product training.

For events, webinars, trainings & individual workshops visit us at www.pedalo.de/academy

We respond to your specific wishes. Let us advise you and inform you about upcoming offers at academy@pedalo.de and +49(0) 73819357-0.

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Subscribe to the newsletter now and
save 10% on your next purchase.

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